MEDSan® COVID-19 IgM/IgG Rapid Test Hybrid

INTENDED USE

MEDSan® COVID-19 IgM/IgG Rapid Test Hybrid is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgM and IgG antibodies against novel coronavirus SARS-CoV-2 infected cells in human whole blood, serum or plasma. This test provides only a preliminary test result. Therefore, any reactive specimen with the MEDSan® COVID-19 IgM/IgG Rapid Test Hybrid must be confirmed with alternative testing method(s) and clinical findings.

INTRODUCTION

Coronaviruses are enveloped RNA viruses that are distributed broadly among humans and that cause respiratory, enteric, hepatic, and neurologic diseases. Seven coronavirus species are known to cause human disease. Four viruses - 229E, OC43, NL63, and HKU1 - are prevalent and typically cause common cold symptoms in immunocompetent individuals. The three other strains - severe acute respiratory syndrome coronavirus (SARS-CoV), Middle East respiratory syndrome coronavirus (MERS-CoV) and 2019 Novel Coronavirus (COVID-19) - are zoonotic in origin and have been linked to sometimes fatal illness. IgM antibody to 2019 Novel Coronavirus can be detected within 7 days. IgG antibody can be detected within 14 days. The combination of both antibodies IgM and IgG can be detected within 14-21 days after infection. IgG remains positive, but the antibody level drops overtime.

PRINCIPLE

The MEDSan® COVID-19 IgM/IgG Rapid Test Hybrid is a lateral flow immunochromatographic assay. The test uses anti-human IgM antibody (test line IgM), anti-human IgG (test line IgG) and goat anti-mouse IgG (control line C) immobilized on a nitrocellulose strip. The burgundy colored conjugate pad contains colloidal gold conjugated to recombinant COVID-19 antigens conjugated with colloidal gold (COVID-19 conjugates) and mouse IgG-gold conjugates. When a specimen followed by assay buffer is added to the sample well, IgM &/or IgG antibodies if present, will bind to COVID-19 conjugates making antigen antibodies complex. This complex migrates through nitrocellulose membrane by capillary action. When the complex meets the line of the corresponding immobilized antibody (anti-human IgM &/or anti-human IgG) the complex is trapped forming a burgundy colored band which confirms a reactive test result. Absence of a colored band in the test region indicates a non-reactive test result. The test contains an internal control (C band) which should exhibit a burgundy colored band of the immunocomplex goat anti mouse IgG/mouse IgG-gold conjugate regardless of the color development on any of the test bands. Otherwise, the test result is invalid and the specimen must be retested with another device.

MATERIALS SUPPLIED

Test kit | Desiccant | 10 µL Pipette | Safety lancet(s) | Alcohol prep pad(s) | Buffer | Package Insert

MATERIAL REQUIRED BUT NOT PROVIDED

- 1. Centrifuge and Pipette (for plasma/serum only)
- 2. Timer

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30 °C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

WARNINGS AND PRECAUTIONS

- 1. The test kit is for professional use only. Do not use after expiration date.
- 2. This package insert must be read completely before performing the test. Failure to follow the insert gives inaccurate test results.
- 3. Do not use if the pouch is damaged or broken.
- 4. Test is for single use only. Do not re-use under any circumstances.
- 5. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- 6. We suggest wearing protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- 7. Humidity and temperature can adversely affect results.
- 8. Do not perform the test in a room with strong air flow, i.e. electric fan or strong air-conditioning.

SPECIMEN COLLECTION

MEDSan® COVID-19 IgM/IgG Rapid Test Hybrid can be performed using either whole blood, serum or plasma.


- 1. Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.
- 2. Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8 °C for up to 3 days. For long term storage, specimens should be kept below -20 °C. Whole blood collected by venipuncture should be stored at 2-8 °C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens.
- 3. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- 4. If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

MEDSan® GmbH | Brandschöfer Deich 10 | 20539 Hamburg | Germany
Tel: +49 40 38 02 38 38 | Fax: +49 40 228 58 56 99 | Email: info@MEDSan.eu

TEST PROCEDURE

Allow test kit, specimen, buffer and/or controls to equilibrate to room temperature (15-30 °C) prior to testing.

- 1. Remove the test cassette from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
- 2. Place the test kit on a clean and level surface.
- 3. If whole blood is collected by fingerstick, it should be tested immediately. Clean your hands before taking blood and wipe off the first drop of blood, use the second drop of blood for testing.
- 4. Add 10 µL of whole blood/serum/plasma to the specimen well (S) of the test cassette, then add 2 drops (about 60 µL) of sample buffer to the specimen well (S) immediately.
- 5. Wait for the colored line(s) to appear. The result should be read at 10 minutes. Positive results may be visible as soon as 2 minutes. Do not interpret the result after 15 minutes.

INTERPRETATION OF RESULTS

NEGATIVE: If only the C band is present, the absence of any burgundy color in the both T bands (IgG and IgM) indicates that no anti-COVID-19 antibodies are detected in the specimen. The result is negative.
IgM POSITIVE: In addition to the presence of C band, if only IgM band is developed, the test indicates for the presence of IgM anti-COVID-19 in the specimen. The result is IgM anti-COVID-19 positive.
IgG POSITIVE: In addition to the presence of C band, if only IgG band is developed, the test indicates for the presence of IgG anti-COVID-19 in the specimen. The result is IgG anti-COVID-19 positive.
IgG and IgM POSITIVE: In addition to the presence of C band, both IgG and IgM bands are developed, the test indicates for the presence of both IgG and IgM anti-COVID-19 in the specimen. The result is IgG and IgM anti-COVID-19 positive.
INVALID: Control line fails to appear. Incorrect specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test kit. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

LIMITATIONS

- 1. Use fresh samples whenever possible. Frozen and thawed samples (especially repeatedly) contain particles that can block the membrane. This slows the flow of reagents and can lead to high background color, making the interpretation of results difficult.
- 2. Optimal assay performance requires strict adherence to the assay procedure described in this insert sheet. Deviations may lead to aberrant results.
- 3. A negative result for an individual subject indicates absence of detectable anti-COVID-19 antibody. However, a negative test result does not preclude the possibility of exposure to or infection with COVID-19.
- 4. A negative result can occur if the quantity of the anti-COVID-19 antibodies present in the specimen is below the detection limits of the assay, or the antibodies that are detected are not present during the stage of disease in which a sample is collected.
- 5. Some specimens containing unusually high titer of heterophile antibodies or rheumatoid factor may affect expected results.
- 6. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Relative Sensitivity (IgG): 98.0%; Relative Sensitivity (IgM): 97.0%
Relative Specificity (IgG): 99.7%; Relative Specificity (IgM): > 99.9%
The test kit does not interfere with other antibodies caused by different viruses such as Influenza A&B, RSV, ADV, HRV, HIV, HBV, HCV or other human coronaviruses (229E, NL63, OC43, HKU1).

REFERENCES

- 1. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011; 81: 85-164.
- 2. Masters PS, Perlman S. Coronaviridae. In: Knipe DM, Howley PM, eds. Fields virology, 6th ed. Lippincott Williams & Wilkins, 2013: 825-58.
- 3. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016; 24: 490-502.
- 4. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019; 17: 181-192.

INDEX OF SYMBOLS

	See instruction for use		Expiry date
	For in vitro diagnostic use only		Batch number
	Store between 2-30 °C		Manufacturer
	Tests per kit		Keep dry
	Catalog number		Do not reuse
	Keep away from sunlight		

MEDSan® GmbH | Brandschöfer Deich 10 | 20539 Hamburg | Germany
Tel: +49 40 38 02 38 38 | Fax: +49 40 228 58 56 99 | Email: info@MEDSan.eu

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MEDSan® GmbH | Brandshofer Deich 10 | 20539 Hamburg | Alemanha